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In the X-ray fluorescence analysis, crucibles and shells made of platinum and its alloys are used on a large scale. However, their almost ideal physical and chemical properties (especially their oxidation and temperature resistance) do not protect them from a whole number of "platinum poisons", which ultimately can destroy the devices. These "poisons" are usually found in the samples, but also in contact materials with which the Pt devices come into contact. The following are some tips on handling to minimize damage and thus increase the service life of the laboratory equipment.

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The physical and chemical properties of platinum and the platinum group metals (Pd, Ir, Rh) characterize them for use in the chemical laboratory:

- resistant to acids and oxides even at high temperatures, by the inert behavior of platinum are also at temperatures above 1000 °C samples not contaminated,
- high melting range,
- low vapor pressure.

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However, platinum forms alloys (phases) with most metals already at medium temperatures, the melting point of which is substantially lower than that of pure platinum.

These phases (also referred to as "platin corrosion") can form due to external contamination or after contact with these elements and lead to various damage mechanisms. For example, platinum reacts with arsenic to form As₂Pt, which has a melting point of 1500 °C. This compound can form a eutectic composition of 72 at.% Pt, which melts at 597 °C.

Already low concentrations of these „platinum-poisons“ (see below) lead to local temperature sinks, which destroy the temperature stability of the platinum device.

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However, other "impurities" are wanted because they have a positive effect on the properties of the platinum. For example, platinum is a relatively soft material and loses its dimensional stability at high temperatures (from 700 °C.). This can be improved by the addition of rhodium, gold or iridium and the lifetime can be increased. Thus, for example, the positive properties of platinum are retained by the addition of iridium, while temperature stability increases. However, due to the evaporation of the iridium, weight loss must be expected.

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In order to prevent / minimize the unwanted contacts with "platinum poisons", different recommendations for the handling result. However, the general rule is that the basic cleanliness of the workplace and the tools and equipment used must be ensured. For example, it is comprehensible that the coatings on which the plating devices are annealed (triangles, grids, etc.) are to be kept free of impurities. These pads are usually coated with platinum, made of platinum, or are covered with oxide-ceramic tubes.

Furthermore, no fluids may enter into the platinum shoes on tools (tongs, tweezers), which would subsequently lead to corrosion. The contact with iron compounds is generally to be avoided.

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The contamination mentioned leads to the diffusion of low-melting elements at high temperatures, whose temperature behavior can lead to destabilization of the device. Accordingly, the use of tongs and tweezers, the tips of which are protected by platinum shoes, should be self-evident.

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Low-melting metals (as mentioned above Arsen), like tin, lead, silicon, antimony, arsenic, boron, phosphorus or bismuth, form with platinum (local) alloys. As a result, in these local areas, the melting point can be exceeded and the device destroyed, or on the other hand cracking by embrittlement at the grain boundaries of these elements. Especially when these organic substances contains phosphorous (e.g., flour).

Phosphor 15 P 30,973761 44 2,19 280	Bor 5 B 10,812 2079 2,04 2550
---	---

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The use of silite bar chamber furnaces for the ash processing involves the risk of contact of the plating devices with exposed silicon carbide heating rods or the contact with their silicon deposits, for example on the base plate. The silicon subsequently diffuses into the platinum lab ware and destroys them by silicon corrosion. At temperatures below 1.000 °C, Pt-Si compounds can be melted to destroy the crucible. Therefore never place Pt-devices on a SiC-base, but to use bases from aluminum oxide or quartz-material. A further risk arises from the SiC heating elements, whose carbide particles can adhere to the crucible wall or contaminate the furnace floor and lead to the same problem as mentioned above.

Conclusion: Keep sufficient distance between heating elements and crucible, clean the oven floor before each sample.

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As is known, the chemical compounds of elements can be dissolved under reduced conditions and thermal influence, and thus harmful elements can be released. In order to avoid a reduction of samples, which must be annealed or incinerated in an oxidizing atmosphere, the preparations should be carried out in open crucibles.

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
The mentioned reduction of samples during thermal processes also takes place by the presence of organic substances, in particular carbon. It is also known that carbon diffuses into the material at the grain boundaries of the platinum, leads to their separation and thus creates porosities (Fig.2).

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Fig.2. At the Pt grain boundaries diffused carbon on the platinum wire

Source: D.Lupton, J.Merker, F.Schölz, On the Correct Use of Pt in RFA-Laboratories, 5th Heraeus User Meeting 3/1990

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Conclusion: It is recommended, for example

- set the flame of the Bunsen burner or other disintegrators as far as possible oxygen-rich, to thereby prevent the release of elements which damage platinum,
- avoid any contact with carbon (deposit not on graphite plates)
- the lowest temperature of the Pt crucible should be used if the contact is indispensable.

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Although phosphorus with platinum forms a eutectic temperature of only 588 ° C, phosphorus damage is generally limited. This is due to the fact that phosphorus rarely occurs in pure form. Figure 3 shows the usual appearance of surface roughening caused by phosphates.

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Fig.3. Formation of brittle platinum phosphide at the platinum-grain boundaries after a phosphate glass melt

Source: D.Lupton, J.Merker, F.Schölz, On the Correct Use of Pt in RFA-Laboratories, 5th Heraeus User Meeting 3/1990

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instructions for handling with platinum labware

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To minimize the damage caused by phosphates, the following information should be implemented:

- *strongly oxidizing conditions during the melt*
- *use of the smallest possible melting temperature*
- *mix the melt sample with the digestant before the melt to reduce the sample concentration prior to contact with the platinum.*

If there is the possibility of carbonizing the samples at a maximum temperature of 900 ° C (organic samples), AuPt10 shells can be used due to their low phosphorus corrosion. However, a lower mechanical strength must also be accepted.

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Apart from platinum, iridium and gold, sulfur reacts very rapidly with non-metals at elevated temperatures. Thus, sulfides in particular lead to brittle dying surfaces, as may occur, for example, in the digestion of glassy tablets for RFA when the sample has a high sulfide content.

Plates can, as a rule, exceed several hundred samples without damage. The contact with sulfur, however, considerably reduces the service life, as the following example shows:

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Fig.5 PtAu5 crucible and cracks due to Pt-S eutectic

Source: D.Lupton, J.Merker, F.Schölz, On the Correct Use of Pt in RFA-Laboratories, 5th Heraeus User Meeting 3/1990

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Crucibles 1 and 2 were used to digest blast furnace slag with 1-2% sulfur content, crucible No. 1 with 27 samples (without cracks) and crucible No. 2 with 50 samples. For comparison, crucible No. 3 with 650 samples used for cement digestions in which no sulfur was present.

Conclusion:

Since the Pt-S eutectic has a very high melting point of 1,240 ° C., damage is relatively straightforward. It already helps if samples with suspected sulfur content are digested in separate crucibles to limit the risks to the selected devices.

Danger:

Sulfur can have a particularly problematic effect on PtRh alloys, since rhodium can form a eutectic with sulfur that melts at 925 ° C.

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The decomposition under oxidizing conditions is particularly important in samples with bismuth, since bismuth oxides are hardly thermodynamically stable. Otherwise, there is a risk of Pt-Bi eutectics, which can lead to the destruction of the crucible with a melting point of 730 ° C. as can be seen in the following figure

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Fig.6 Cut through a Pt-Bi eutectic in a Pt-crucible wall
Source: D.Lupton, J.Merker, F.Schölz, On the Correct Use of Pt in RFA-Laboratories, 5th Heraeus User Meeting 3/1990

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Copper is in itself not a major hazard to plating devices, if attention is paid to oxidizing conditions during digestion. Otherwise, connections can occur at the grain boundaries of the platinum, which lead to embrittlement of the surface or even breakage.

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Further damage is caused by alkaline samples, acids, salts and halogen compounds. Examples:

- *even at room temperature, platinum is chemically attacked by strongly oxidizing acid mixtures (aqua regia, halogen compounds)*
- *since alkalies (A. metal hydroxides and cyanides) act as oxygen carriers, platinum can oxidize to yellow to brown platinum oxide during these digestions, potassium compounds reacting more intensively than sodium compounds (concerns, for example, soda and soda potash)*
- *In the long term, platinum forms an oxidation layer at room temperature, which evaporates at high temperatures. Partially, on the other hand, the alloying of rhodium helps.*

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It should also be noted that after the absorption of hydrogen by platinum at 400 ° C, this diffuses through the device wall at higher temperatures and thus reduces the annealing sample.

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Intensive cleaning can be carried out by melting with potassium pyrophosphate or potassium hydrogen sulfate, if boiling with suitable solvents does not lead to the desired result.

Melting residues can be abraded with abrasives (corundum) (do not use carbonaceous or metal carbide-containing materials), which must be removed afterwards with hydrofluoric acid.

Metal layers on circuit boards are peeled off with analytic acids, cleaned with deionized water and stored dry. For this a desiccator is recommended.

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In summary, the following information to increase the useful life of the platinum lab ware (without claim to completeness):

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Quelle: Heraeus

- *never touch platinum instruments with a bare hand*
- *do not place Pt lab ware on SiC basis*
- *keep distance from SiC heating elements*
- *place Pt crucible only in the burning flame to avoid contact with unburned gas, ensure oxygen-rich flame adjustment*
- *use Pt shoes for instruments / tools, especially for hot appliances*
- *use separate crucibles for different (aggressive) samples*
- *ensure oxygen access (oxidizing conditions) for carbonaceous / organic samples (add small amounts of ammonium nitrate before heating)*
- *avoid mechanical damage to the crucibles*
- *always treat alkaline samples in the covered crucible, whereby the carbon dioxide released during digestion acts as protective gas and protects the crucible from oxidation*
- *annealing as a cleaning method damages the Pt lab ware by diffusion of impurities into the platinum*
- *in the case of simultaneous annealing of several crucibles in the muffle furnace, they can weld each other when they are in contact*
- *dispersion-enhanced (DV) plating devices can have an extended lifetime by a factor of 3, particularly at high temperatures*
- *for unknown samples already used very used crucibles*

The sample-dependent use of crucibles also has significance for the selection of the alloy. Here are some hints for the selection of the crucible alloy as a function of the sample material:

- *Pt for simple routine digestions*
- *PtAu5 for low wetting*
- *AuPt10 for phosphorous-containing samples, but their temperature is limited to about 900 ° C*
- *never use Rh-containing alloys for sulfur-containing samples.*

For further information we are at your disposal.

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- H.Jehn: "High Temperature Behaviour of Pt-Group Metals in Oxidizing Atmospheres", Journal of the Less Common Metals, 100 /1984), 321

Source:

- J.Merker, F.Scholz, D.F.Lupton: "On the correct use of platinum in RFA laboratories", 5. User meeting X-ray fluorescence and spark emission spectroscopy, Dortmund, 1998

- H.Jehn: "High Temperature Behavior of Pt-Group Metals in Oxidizing Atmospheres", Journal of the Less Common Metals, 100/1984), 321

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The process of recycling begins with the expert advice on the effective treatment of your waste in order to evaluate the costs and the yield. Precious metals are found in many different products across the industry, for example in fractured gold, old jewellery, scratches, crucibles, targets, wires, contacts, electroplating waste, electronic scrap, shields and diaphragms, punching and casting waste, filing and cutlery, dental waste (crowns and bridges) As a result, the assignment to a processing technology already determines the effectiveness of the recycling process. You support us in this decision by providing information about your material. Please use our „Checklist for recycling material“.

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After deciding on the kind of preparation of the material, it runs individually per customer and batches all sections of the work-up, from the homogenization (melt, preparation, analysis) to the provision of the fine metal account. The precious metal can be provided virtually on a metal account (see section "The metal account as a central part of your precious metal management") for the next demand, physically delivered in the form of bars or granules or paid out at the appropriate exchange rate.

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Our service for you:

- „transparent melt“ - you are present during melting and analyzing your material on request
- precious metal management including collection of the old material
- sampling and efficiency calculation
- individual billing
- free metal account management and transfers within Germany

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Our strength lies in customer service. We take care of our customers personally, deal with problems, grow together with tasks and always find specific solutions for your needs.

We would be happy if we could prove it. Just send us your non-binding inquiry and rate it yourself.

Note: You can also use our website for information and orders: www.mk-edelmetall.de

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Telefon: +49 (0) 36424-8110
E-Mail: mail@mk-edelmetall.de
www.mk-edelmetall.de