

Laborgeräte für die Röntgen-Fluoreszenz-Analyse(RFA)
labware for x-ray fluorescence analysis (XRF)

Tiegel für Katanax
crucibles for Katanax





Tiegel für Katanax

crucibles for Katanax

Material z.B.: Pt/Au 95/5; DV Pt/Au 95/5

Bezeichnung description	ø innen internal øU [mm] øO		Höhe height H [mm]	ø Gesamt total øA	Bodendicke base thickness B [mm]	ca. Gewicht approx. weight [g]	m&k-code
ohne Randverstärkung without reinforced rim	28,0	40,0	32,0	40,5	0,4	30,0	TKA4
ohne Randverstärkung without reinforced rim	28,0	39,5	32,0	41,0	1,0	75,0	TKA5
mit Randverstärkung with reinforced rim	28,0	39,5	32,0	44,0	0,4	30,0	TKA6



mit Randverstärkung
with reinforced rim



ohne Randverstärkung
without reinforced rim

Abgießschalen für Katanax

casting dishes for Katanax

Material z.B.: Pt/Au 95/5; DV Pt/Au 95/5

ø Schalenrand flange dish K [mm]	ø innen internal øU [mm] øO		Höhe height H [mm]	Bodendicke base thickness B [mm]	ca. Gewicht approx. weight [g]	m&k-code
49,0	40,0	42,0	7,0	1,0	48,0	ASKA17
49,0	40,0	42,0	4,5	1,0	45,0	ASKA18
41,0	32,0	33,0	6,0	1,0	35,0	ASKA19

