

Laborgeräte für die Röntgen-Fluoreszenz-Analyse(RFA)
labware for x-ray fluorescence analysis (XRF)


Tiegel mit Doppelkonus
crucible with double-cone



Tiegel mit Doppelkonus
crucible with double-cone

Material z.B.: Pt/Au 95/5; DV Pt/Au 95/5


σ zwischen Konen <i>σ between two cones</i> [mm]	σ innen <i>internal</i> σU [mm] σO	Höhe <i>height</i> H [mm]	Bodendicke <i>base thickness</i> B [mm]	ca. Gewicht <i>approx. weight</i> [g]	m&k-code
34,0	29,0 38,0	39,3	0,8	53,0	TDK1
36,0	34,0 40,0	40,0	0,8	64,0	TDK2
41,0	36,0 45,0	43,3	0,8	75,0	TDK3



Die Perl-Höhe (Ph in der Skizze) beträgt standardmäßig 3 mm. Andere Perl-Höhen können auf Wunsch angefertigt werden.
 The Perl-height (Ph in the sketch) is 3 mm as standard. Other Perl-heights can be made on request.

Tiegel mit drei Zapfen
crucible with three Pins

σ innen <i>internal</i> σU [mm] σO	Höhe <i>height</i> H [mm]	Bodendicke <i>base thickness</i> B [mm]	ca. Gewicht <i>approx. weight</i> [g]	m&k-code
40,0 30,0	30,0	0,3	20,0	TDZ1
40,0 30,0	30,0	0,45	33,0	TDZ2
44,0 32,0	30,0	1,0	60,0	TDZ3



m&k gmbh
 Im Camisch 49
 07768 Kahla

Telefon: +49 (0) 36424-8110
 E-Mail: mail@mk-edelmetall.de
www.mk-edelmetall.de